


1. TLC method development




Mobile phase:
75% Hexane / Ethyl Acetate 25%

Compound	Rf	CV
1	0.44	2.26
2	0.34	2.94

$\Delta CV = 0.68$

2. Choice of the column according to the ΔCV & crude sample mass

Crude sample: 32mg
Column: PF-15SIHC-F0004
Loading capacity: 0.8%



Customer has chosen to use PF-15SIHC-F0004 column to obtain a better separation (efficiency & purity) than with a IR-50SI-F0004 column.

3. Flash conditions

Device: puriFlash® 4125-iELSD (or now puriFlash® 5.125 pack-iELSD)

Solvents: A: Hexane

B: Ethyl Acetate

Column: PF-15SIHC-F0004

Flow rate: 7mL/min

Injection mode: Solid deposit with celite (Dry-load F0004)


Crude sample: 32mg

Detection: UV Scan 254-280nm (orange), iELSD (T°: 35°C, Gain: 6) (purple)

Mode: Automatic Gradient Optimization

Elution conditions:

CV	A (%)	B (%)
0	100	0
2	100	0
12	85	15
36	85	15



To achieve this purification:

You will need

- puriFlash® 5.125
[Discover it](#) [Add to card](#)
- Integrated ELSD
[Discover it](#) [Add to card](#)
- puriFlash® column PF-15SIHC-F0004
[Discover it](#) [Add to card](#)
- puriFlash® Dry-load PF-DLE-F0004
[Discover it](#) [Add to card](#)

We highly recommend

- Trolley
AYHF20 [Add to card](#)
- Manometer ELSD
FJ6720 [Add to card](#)
- Tube holding claw 16mm
AYHECO [Add to card](#)

Download our App

"TLC to Flash & Prep Chromatography"
to make your TLC developments easier
and faster.

