


1. TLC method development




Mobile phase:
50% Heptane / Ethyl Acetate 50%

Compound	Rf	CV
1	0.66	1.52
2	0.43	2.33
3	0.26	3.8

$\Delta CV_{2-1} = 0.81$

2. Choice of the column according to the ΔCV & crude sample mass

Crude sample: 500mg
Column: PF-15SIHP-F0025
Loading capacity: 2%



Customer has chosen to use PF-15SIHP-F0025 column to obtain a better separation (efficiency & purity) than with a IR-50SI-F0040 column.

3. Flash conditions

Device: puriFlash® 4250-iELSD-MS
(or now puriFlash® 5.250-MS pack iELSD)

Solvents: A: Heptane
B: Ethyl Acetate

Column: PF-15SIHP-F0025

Flow rate: 20mL/min



Injection mode: Liquid injection

Injection volume: 1 mL

Detection: UV 254nm (Black), ELSD (Orange),
MS (APCI): TIC, XIC (m/z 244) (Blue)

Eluent conditions:

t (min)	A (%)	B (%)
00:00	95	5
02:00	95	5
25:00	10	90



To achieve this purification:

You will need

- puriFlash® 5.250
[Discover it](#) [Add to card](#)
- puriFlash® MS
[Discover it](#) [Add to card](#)
- Integrated ELSD
[Discover it](#) [Add to card](#)
- puriFlash® column PF-15SIHP-F0025
[Discover it](#) [Add to card](#)

We highly recommend

- Trolley
AYHF20 [Add to card](#)
- Manometer ELSD
FJ6720 [Add to card](#)
- Safety waste cap with container 5L + Filter
B1SUJ0 [Add to card](#)

Download our App

"TLC to Flash & Prep Chromatography"
to make your TLC developments easier
and faster.

