

1. Analytical conditions

Solvents: A: Water

B: Acetonitrile (0.1% Formic acid)

HPLC Column: Waters Acquity BEH C18 250x4.6mm 1.7 μ m

Flow rate: 1mL/min


Injection mode: Liquid injection

Injection volume: 10 μ L

Concentration: 0.5mg/mL in Acetonitrile

Detection: UV 315nm

Elution conditions: 40% Water / Acetonitrile 60%
(Isocratic mode)



2. Prep conditions

Device: puriFlash® 4250 (or now puriFlash® 5.250)

Solvents: A: Water

B: Acetonitrile (0.1% Formic acid)

Column: PF15C18HP-150/212

Flow rate: 20mL/min

Injection mode: Liquid injection

Crude sample: 100mg


Detection: UV 315nm (black),

Scan 200-600nm (orange)

Pressure: 7bar

Elution conditions:

t (min)	A (%)	B (%)
00:00	40	60
23:00	40	60
24:00	0	100
29:00	0	100



To achieve this purification:

You will need

- puriFlash® 5.250
[Discover it](#) [Add to card](#)
- puriFlash® column PF15C18HP-150/212
[Discover it](#) [Add to card](#)

We highly recommend

- 1 mL Stainless Steel Loop with RFID tag
AYHET0 [Add to card](#)
- Ballasting for 1/8" tubing - 5 units
DZ7360 [Add to card](#)
- Safety solvent caps kit - 4 units
BODANO [Add to card](#)

