



puriFlash® R-Vap

Your laboratory rotary evaporator

Super intuitive. Ultra-efficient.
Simply indispensable

Components

Vacuum rotary evaporator with evaporation flask, receiver flask and cooler.

Features:

- Rotation speed range 20-300 RPM
- Rotation: Clockwise / Counterclockwise
- Lift mode: Manual / motorized
- Immersion angle: 0° to 60°
- Motor power: 40W


Heating Water Bath.

Features:

- Water (ambient temperature up to 100°C) - Oil (ambient temperature up to 230°C).
- Volume Capacity: 5L
- Temperature control mode: PID Digital display
- Timer: yes
- Overheat cut-off protection temperature: yes
- Heating power: 1300W

Glassware dimensions (standard version or plastic coated)

- Evaporation flask: 50mL, 100mL, 250mL, 500mL, 1L, 2L, 3L.
- Receiving flask: 50mL, 100mL, 250mL, 500mL, 1L, 2L.
- Vertical condenser column: 1200cm², 1500cm².



Safety

Liquid bath safety shield for puriFlash® R-Vap



Power Supply

Voltage: 110V ; 220V
Choice of power cord: European, UK, US, Switzerland, Japan, India, China

Dimensions & Weight

Weight: 19kg

